AFFINIS® DCode heavy body System 360

DE	Gebrauchsinformation	3	PL	Instrukcja użycia	55
EN	Instructions for use	7	SL	Navodila za uporabo	59
FR	Mode d'emploi	11	LT	Naudojimo instrukcijos	63
ES	Información para el uso	15	ET	Kasutusjuhend	67
IT	Informazioni per l'uso	19	CS	Návod k použití	71
SV	Bruksanvisning	23	HU	Használati utasítás	75
NL	Gebruiksaanwijzing	27	SK	Návod na použitie	79
DA	Brugsanvisning	31	BG	Инструкции за употреба	83
NO	Bruksanvisning	35	LV	Lietošanas instrukcijao	87
FI	Käyttöohje	39	TR	Kullanma talimatı	91
EL	Πληροφορίες χρήσης	43	RU	Инструкции по применению	95
PT	Instruções de Utilização	47	ZH	使用说明书	99
RO	Instrucțiuni de utilizare	51			

Definition

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EN

Definition

AFFINIS DCode System 360 heavy body is a silicone-based scannable dental impression material indicated as a tray material for use with commercially-available automated mixing machines.

Material type

Polyvinylsiloxane, addition-type, surface activated silicone elastomer.

heavy body System 360

Base: yellow
Catalyst: caramel
– ISO 4823, Type 1, high consistency

Indications

- One step putty/wash impressions (double mix)
- Two step putty/wash impressions
- Functional impressions
- Transfer coping impressions / pick-up impressions
- Relining impressions
- Impressions for implantology (after healing is complete)
- Digitalisation of the impression with the aid of a dental scanner without additional surface treatment (only in combination with AFFINIS DCode materials)

Contra-indications

Provided that the product is used as directed, there are no known contra-indications.

Side effects / Interactions

Polyvinylsiloxanes have a very good biological compatibility and up to now no harmful reactions or secondary effects on patients and/or

dental personnel are known. Impression materials are intended to set in the mouth of the patient, however, they should remain in the mouth not more than twice the recommended setting time. Although they have reasonably high tear strength, care should be taken that no portion of the impression material remains in the interdental spaces or in the sulcus. Undercuts should, in certain instances, be blocked out before taking the impression.

Trays

Tray selection depends on impression technique and professional preference of stock trays or custom trays. For perfect adhesion, we recommend applying a thin layer of COLTENE Adhesive or any other brand of adhesive specified for use with polyvinylsiloxane impression materials.

Preparation

1. Hold the SYS360 Cartridge facing upwards with the outlet channel closed.
2. Gently clip the safety cap of the outlet channel while holding it upwards; then remove and discard the safety cap.
3. Insert the SYS360 Cartridge according to the manufacturer's instructions for use of the mixing machine.
4. Start the mixing machine and extrude material out of the outlet channel without the mixing tip until equal amounts of base and catalyst are dispensed. This ensures that the plungers are at the same height to achieve an optimal mix.
5. Carefully wipe off extruded material using a paper towel to avoid contamination between the base and catalyst in the outlet channels.
6. Attach a yellow dynamic mixing tip onto the outlet channel of the SYS360 Cartridge. Pay attention that the hexagon section of the

yellow dynamic mixing tip is correctly aligned onto the drive shaft.

7. When the yellow dynamic mixing tip is correctly placed, the yellow fixation ring must be positioned over the dynamic mixing tip and turned right (clockwise ¼ turn) until it is securely fastened. The mixing components inside of the yellow dynamic mixing tip must rotate when the machine is started.
8. Before every use reject the first 3 cm of the mixed impression material and press out the following material without interruption.
9. Fill the mixed impression material directly into the selected impression tray. In order to prevent the formation of air voids, the yellow dynamic mixing tip should always remain immersed in the material.
10. In order to prevent excess material coming out, the plungers of the mixing machine should be disengaged after extruding material.
11. The yellow dynamic mixing tip should be left on the SYS360 Cartridge until the next application, since it acts as a seal. It should be wiped with a disinfection solution.
12. To remove the used yellow dynamic mixing tip, the fixation ring must be turned counter-clockwise and removed. The fixation ring is reusable.
13. Before the next application, check to ensure that the outlet channels of the SYS360 Cartridge are not clogged and that equal amounts of material can be extruded. Proceed according to point 5 as follows.
14. Remove the empty SYS360 Cartridge out of the mixing machine and discard it according to local regulation.

Notice

SYS360-Cartridge

- Do not drop the SYS360 Cartridge, since this may damage it, leaving it defective.

Yellow dynamic mixing tip

- The yellow dynamic mixing tip should be left on the SYS360 Cartridge until the next application, since it acts as a seal.
- To guarantee trouble-free application, we recommend that the Coltène/Whaledent yellow dynamic mixing tips be used exclusively.

Recommendations for use

One step putty/wash technique (double mix)

While loading the tray, the dentist may begin to apply the wash-material around the preparation. Tray loading and application of wash-material around the preparation must be completed at the same time. Keep the oral tip immersed in the material at all times to avoid trapping air. *Immediately place the filled tray into the patient's mouth. Press for 2-3 s and hold in position until the material is completely set.*

"Two step technique"

If a 2-Step «heavy body/wash» technique is used, the primary impression made with heavy body or putty-material, must be carefully cleaned and dried before further use. This will help guarantee good adhesion between the wash- and tray-material.

For extended working time, material can be placed in the refrigerator prior to use. High temperatures will speed up the setting process and low temperatures will slow it down.

Important

Always wear gloves.

Skin secretions, latex gloves and surfaces contaminated by latex gloves can interfere with the setting process of polyvinylsiloxanes. The material and surfaces where the impression will take place (teeth, preparations, retraction cords, etc.) should only come into contact

with thoroughly washed and rinsed gloves (wash 15 s with a detergent, rinse with lukewarm tap water for another 15 s) or vinyl-/nitrile gloves. Compounds containing eugenol or hemostatic substances can also impede the setting process. If (H₂O₂) hydrogen peroxide is used for disinfection, it is recommended to thoroughly rinse with lukewarm water to avoid the formation of air voids.

Check intraorally if the material is completely set before removing from the mouth.

Disinfection

The impression should be rinsed under running tap water after removal from the mouth. After rinsing, disinfection with a suitable commercial dental disinfectant solution (according to the manufacturer's recommendations) which will not affect the impression surface or dimensions. Acrylic trays should be protected against water absorption.

Optional

The following points should be observed when autoclaving impressions:

1. Only use autoclavable components (e.g. PRESIDENT Tray AC).
2. Rinse and clean impression thoroughly under flowing, lukewarm water.
3. Impressions can be autoclaved directly after taking the impression.
4. Autoclave in steam steriliser at 134°C/273°F (Prion program)

Important

When sterilising implant impressions, be sure clarify in advance with the manufacturer, whether or not the implant components (e.g. impression posts, etc.) are autoclavable.

Scanning of the impression

The material can be digitalised with conventional dental scanners without the additional use of scan powders or sprays. The quality of the representation depends on the impression scanner used. The recommendations of the respective manufacturers on preparation of the teeth and the procedural instructions for optical recording are to be followed.

Fabrication of models

The stone models can be poured after 30 min at the earliest. The impression remains dimensionally stable for a practically unlimited period of time (for at least 1 week). The surface tension will be reduced and pouring will be facilitated if the impression is briefly washed out with a detergent and rinsed thoroughly in lukewarm clear water afterwards. COLTENE impression materials are compatible with all dental stones, epoxy resin and polyurethane.

Electro plating

The impression is compatible with all commercially available copper- and silver plating baths.

Cleaning of trays

The impression can be removed with a blunt instrument. Soaking in a universal commercial solvent or petroleum ether will dissolve the COLTENE Adhesive. These solvents should only be used in a well-ventilated area. The trays can also be ultrasonically cleaned and sterilized.

Shelf life and storage

AFFINIS DCode fulfils the intended purpose at least until the expiry date, which is displayed on the containers, in well sealed containers at temperatures 15–23 °C / 59–73 °F and 50 % relative air humidity. Im-

pressions should be stored at normal room temperature, avoid exposure to heat and sun.

Caution

Federal law restricts this device to sale by or on the order of a dentist.

Marking

The expiry date and **LOT**-number are displayed on the package and cartridge.

Technical data ISO 4823

Measurements are based on 23 °C / 73 °F room temperature and 50% relative humidity.

AFFINIS DCode heavy body System 360

Mixing time (15ml):	0:15 min
Total working time:	0:45 min
Time in mouth (35 °C / 95 °F):	2:00 min

Date of issue

01-2017

Only supplied to dentists and dental laboratories or upon their instructions.

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Coltène/Whaledent AG